



PRE-POST 1/2%®

A SANITIZING PRE & POST DIP

New Item #	Case UPC Mfg. No. 071206	Each UPC Mfg. No. 071206	Case Size	Case Weight	Cases /Pallet
#1244010	144010	144317	4 x 1 gallon case	37 pounds	48
#1244058	144058	N/A	5 gallon pail	47 pounds	36
#1244157	144157	N/A	15 gallon drum	138 pounds	9
#1244553	144553	N/A	55 gallon drum	493 pounds	4

A sanitizing pre & post iodine teat dip, NPE-free

- ½% Iodine, 5% glycerin
- An aid in reducing the spread of organisms which may cause mastitis
- For use before and after milking
- Emollient system provides excellent skin conditioning properties, minimizes teat irritation
- NPE-Free to reduce environmental impact
- FDA regulated, NDC #060282-2106

Directions

For Pre-Milking

1. Wash each individual teat with a sanitizing udder wash.
- NOTE:** For best results, keep hair around teat base clipped.
2. Hand strip teats to initiate let down.
3. Pre-dip entire teat with Pre-Post 1/2%® on all four quarters. Repeat process on next 4 to 6 cows or until first cow prepped has at least 1 minute contact time.
4. Wipe each individual teat with a clean paper towel. Utilize a “pulling down” motion which removes iodine solution, germs, and sediment from teat skin. Very important!
5. Place machine on cow and repeat wiping process on remaining cows that have been pre-dipped.

For Post Milking

Dip entire teat with undiluted dip solution. Allow to air dry and remain on teats for prolonged contact time. Clean, sanitize, and thoroughly dry each teat prior to next milking. Replace dip solution when visibly cloudy or dirty.

Chemical Characteristics

Appearance	Dark Brown Liquid
Biodegradable	Yes
Corrosive	No
DEA-Free	Yes
Fragrance	Iodine Odor
Non-Flammable	Yes
NPE-Free	Yes
NSF Certified	No
pH (ready-to-use)	4.8-5.2
Phosphorus	No
VOC Compliant.....	None

Ingredients

The active ingredient in Pre-Post ½%® is 0.5% available iodine. Other ingredients include: iodophor complex, 5% glycerin, nonionic surfactant, xanthan gum, thickening agent, buffering agent, and pH adjusters.



STEARNS PACKAGING CORPORATION

4200 Sycamore Avenue (53714) • P.O. Box 3216, Madison, Wisconsin 53704-0216

Telephone: (800) 655-5008 • (608) 246-5150 • Fax: (608) 246-5149

E-mail: stearns@stearnspkg.com • Website: <http://www.stearnspkg.com>